

XTRACTIVE[®]

POLARIZED™

THE ONLY AND BEST
EVER PHOTOCROMIC
POLARIZED LENS¹



QUICK FACTS

Transitions® XTRActive® Polarized™ lenses:

- Expand the *Transitions XTRActive* product line with the additional benefit of polarization
- Start clear indoors, darken in the car, and achieve extra darkness and polarize outdoors in the sun
- Provide extra protection for eyeglass wearers who are frequently in intense light situations or who are very light sensitive
 - The unique polarization technology reduces glare outdoors, providing sharper vision, a larger field of view and bright, vivid colors

NEW TECHNOLOGY

- Use a breakthrough, exclusive, multi-layer matrix which hosts new *Transitions XTRActive* broad-spectrum dyes for more darkness with new ultra-fast dichroic dyes that selectively absorb light waves in one direction for polarization

SPECIFICATIONS

- **Activate behind the windshield**, which is normally a barrier for entry into the category for eye care professionals when recommending
- **Polarization:** Outside, *Transitions XTRActive Polarized* lenses achieve a polarization efficiency of up to 90 percent, which is similar to polarized sunglasses
- **Fadeback Speed:** The lenses fadeback up to 2 times faster than *Transitions® Vantage®* lenses
- **Clarity indoors:** Filters 2.6 times more harmful blue light than clear lenses indoors

¹ Compared to clear to dark photochromic lenses